

# made in Japan made by OSADA

We're one with the doctor  
and the dental hygienist.

**DAISY 2**

 **OSADA**  
Refine for the future

# Anytime, Anywhere, Anyone.

Like the daisy's flower language, "I feel the same way you do," Osada portable unit Daisy is considerate of the patient's feelings and treats them gently, as if they were a smiling flower. It is designed to be patient- and doctor- friendly, and can be easily transported anywhere, anytime, to provide the same level of care as in a clinic room, even when visiting a patient.



---

DAISY 2

---





Osada's home-visit series provides reassuring treatment even in the home, featuring not only functionality but also a gentle, compact form that avoids intimidating patients. We conduct regular visits to patients who have difficulty traveling to the clinic, delivering smooth, efficient treatment with no wasted time. Whether the patient is in a wheelchair or seated in a home chair, care is provided in a posture that minimizes strain and maximizes comfort.



Bring the clinic to the home.



• For storage



• For solo-operation by a dentist



• For dual-operation by a dentist and an assistant

3 installation options

The slim, compact design allows the unit and vacuum to function as a single space-saving unit. They can also be separated when treatment requires two operators (a dentist and an assistant). Both the unit and vacuum can accommodate support columns for adjusting the height of the instrument holder, enabling comfortable treatment while standing.

High Performance & Operability



- Standard specification Instrument
- ☐ G1 Micromotor (G1M)
  - ☐ 200,000 rpm contra angle handpiece (WG-99LT)
  - ☐ Straight handpiece (SHG02)
  - ☐ Multi-syringe F (MS-F)
  - ☐ Multi-purpose ultrasonic instrument with LED light-Enac 11L (SE-11L)



• High-performance Instrument



• Water tank (400ml capacity)



• Drain tank (280ml capacity)

**High Performance** The instrument holder can be stored in the unit body and does not need to be removed, saving time for setup, operation, and cleanup. The standard specification includes a 200,000 rpm contra handpiece and a straight handpiece, enabling high-quality dental treatment even in patients' homes. The water tank and drain tank are easily detachable for simple maintenance.



• Touch control panel



• Power control switch



• Maintenance operation

Operability

All operations can be performed on the main unit's touch control panel. A streamlined power control switch contributes to a simple, easy-to-understand design—even for doctors using it for the first time. Intuitive controls allow easy adjustment of rotation speed, power settings, and mode selection.

For  
easy transport



• On the carrying cart



• With the shoulder belt



• With the full cover

Mobility  
(Optional)

The carrying cart and shoulder strap make transportation easier and reduce the burden of carrying. The full protective cover prevents scratches and stains, while also providing storage space for small accessories needed in mobile dental treatment.



The reassurance that  
“My dentist will come  
to see me.”

With the improved performance of the portable unit, this single unit can transform any environment into a fully functional dental clinic. From an era when patients came to you, to an era where you visit them. This unit provides reassurance as it helps you build a practice system that meets the needs of changing times and reflects a patient-centered approach.





OSADA ELECTRIC CO.,LTD.  
Since 1935

#### OSADA PORTABLE UNIT

The Osada portable unit was first launched in 1984 with a single goal : “ to bring clinic-level treatment out of the clinic room.” Looking ahead, dentists will increasingly visit patients who cannot come to the clinic to provide care. For more than 30 years, we have pursued in-home treatment so that the treatment room can be brought directly to the patient’s home. The evolved portable unit is used not only for home dental visits but also for dental treatment and oral care in disaster-affected areas without compromising functionality.

Designation : Portable/Optional dental unit    Product name(Model)    : Osada portable unit Daisy(OPU-D2)

● Please be sure to read the attached document and instruction manual before using the product.

## Osada Electric Co., Ltd.

5-17-5, Nishi-gotanda, Shinagawa-ku, Tokyo 141-8517, Japan  
<https://osada-group.jp/global/index.html>

[Click here for Osada website](https://osada-group.jp/global/index.html)



[Contact us here](#)



\*The colors may appear different from the actual colors due to the printing ink. \*Please note that specifications are subject to change without notice for improvement.