

# made in Japan made by OSADA

We're one with Dr. & Dental Hygienist

**DAISY 2**

 **OSADA** Since 1935

Any time,  
Any where,  
Any one.

Like the daisy's flowery language, "I feel the same way you do," Osada portable unit Daisy is considerate of the patient's feelings and treats them gently, as if they were a smiling flower. It is designed to be patient- and doctor-friendly, and can be easily transported anywhere, anytime, to provide the same level of care as in a clinic room, even when visiting a patient.



---

DAISY 2

---



Osada's home-visit series offers not only functionality, but also a gentle, compact form that is not intimidating to patients, providing them with peace of mind when they are at home. Regular visits to the homes of patients who have difficulty visiting the dental clinics can be performed smoothly and effectively.

Patients can receive dental treatment in a wheelchair or in a table chair at home, in a comfortable posture that is easy for them to bear.



## High Performance & Operability



### Standard specification Instrument

- Micromotor G1 motor (G1M)
- 200,000 rpm contra angle handpiece (WG-99LT)
- Straight handpiece (SHG02)
- Multi-syringe F (MS-F)

### Option

- Multi-purposed ultrasonic instrument Enac 11L (SE-11L) with LED light.



• High-performance Instrument



• Water tank (400ml capacity)



• Drain tank (1,100 ml capacity)

## Treatment room carried directly to patient's home



• For storage



• For solo-operation by a dentist



• For dual-operation by a dentist and an assistant

### 3 ways for installation

The slim and compact design allows the unit and vacuum to be used as one unit to save space. It can also be separated according to the situation when treating by two with an assistant. Both the unit and the vacuum can be supported by pillars, and the height of the instrument holder can be adjusted so that the treatment can be done in comfortably positioned in a standing posture without strain.

### High performance

The instrument holder can be stored in the unit body and does not need to be removed, taking no time for preparation, operation, and cleanup. 20,000 rpm contra handpiece and straight handpiece which is standard specification, enables high-quality dental treatment even patients' home. The water tank and the drain tank are easily detachable which is designed for easy maintenance.



• Operation touch panel



• Power control switch




• Maintenance operation

### Operability

All operations can be performed on the control touch panel of the main unit. The combination of the power control switch has resulted in a simple design that is easy to understand even for first-time doctors. The intuitive switch allows easy adjustment of rotation speed, power settings, and mode selection.

## High Specification Option




Oral care (Option)

Combination of the vacuum unit and the oral care set of Daisy allows for safe and reliable cleaning of the oral cavity. Warm water and medicament can be used by supplying water from a paper cup. The tip of the oral care brush simultaneously pours water and suction to prevent accidental ingestion.

## For smooth carriage



• On the carrying cart



• With the shoulder belt



• With the full cover


Mobility (Option)

The carrying cart and the shoulder belts can be used for less burdensome carrying. A full cover prevents scratches and stains and also provides storage space for small items required for visit dental treatment.



The reassurance that  
“My dentist will come  
to see me.”

With the improved performance of the portable unit, this one unit can transform any environment into a dental clinic. From the era of having patients come to you to the era of visiting patients. This is a reassuring unit that will help you realize a practice system that meets the needs of the changing times, as well as an attitude of compromise.



OSADA ELECTRIC CO.,LTD.  
Since 1935

#### OSADA PORTABLE UNIT

In 1984, Osada portable unit was launched. "We want to develop a treatment device that can be delivered out of the clinic room." In the future, doctors will visit patients who cannot come to the clinic to provide treatment. For more than 30 years, we have been trying to realize treatment at home so that the treatment room can be carried directly to the patient's home. The evolved portable unit is used not only for visiting dental treatment but also for dental treatment and oral care in disaster-stricken areas without inconvenience.

● Please be sure to read the attached document and instruction manual before using the product.

#### Osada Electric Co., Ltd.

5-17-5, Nishi-gotanda, Shinagawa-ku, Tokyo 141-8517  
<https://osada-group.jp/global/index.html>

[Click here for Osada Group website](https://osada-group.jp/global/index.html)



\*Products may be changed for revision without notice. Contact us or Osada distributors for any questions.  
\*Actual color of products may differ to the printed catalog.